



SomaTherapy-ED



Instant Relief for Erectile Dysfunction.

SomaTherapy ED is designed to produce an erection at the touch of a finger by combining vacuum technology with the natural function of the body.

SomaTherapy ED can be used in conjunction with Erectile Dysfunction medication, or can be used as a supplementary therapy. Please consult your doctor for further information about the range of available treatment options, and to discuss whether SomaTherapy ED is the right treatment for you.

Why choose SomaTherapy ED?

- Unlike medication, SomaTherapy ED is designed to work immediately.
- Patient controlled.
- The SomaTherapy ED package comes with tension bands which are designed to maintain the erection for the duration of intercourse.
- Available in manual or battery operated variations.



Touch II VED System

- Super-quiet battery operation
- Activation at the press of a button
- Ergonomic design for men with limited dexterity or arthritic finger joints

Order code: AUG 15013



Response II VED System

- Designed for one finger operation.
- "Ultra Ease" mechanism ensures the same gentle pressure as battery operation

Order code: AUG 15019



Endotherapeutics

This medical device is approved for use by the TGA.
ARTG Number: 123413

SomaTherapy-ED



Both the Response II and Touch II are portable kits and include the following:


- Clear plastic cylinder.
- Loading cone for application of the tension rings.
- Cylinder inserts to prevent scrotal tissue from being trapped within the cylinder.
- Tension rings.
- Lubricant

Tension rings available in various models and sizes.


Tension rings and lubricant are considered to be consumable items and can be ordered separately.*

How to use SomaTherapy ED

Apply the system over the penis.



Activate the pumping mechanism. This will draw blood into the penis, causing an erection.



Once the penis is fully engorged, the tension ring is placed at the base of the penis, maintaining the erection.



*Disclaimer:

Tension rings have a limited life expectancy and cannot be expected to be used indefinitely. Their useful life will change depending upon the circumstances of their use. The lifetime warranty does not cover this consumable item.



Endotherapeutics

Suite 1, Level 1
51 Rawson Street
Epping NSW 2121
Australia

T: +61 2 9869 2868

F: +61 2 9869 7723

www.endotherapeutics.com.au

info@endotherapeutics.com.au